

CANDOR The ELISA Experts **APPLICATION GUIDE**

For CANDOR Bioscience Products

Recommended Application	The Blocking Solution™ Cat. No.: #110	LowCross-Buffer® Cat. No.: #100	LowCross HRP® Cat. No.: #200	LowCross® HRP-Stab Cat. No.: #270	Liquid Plate Sealer® Cat. No.: #160	Antibody Stabilizer Cat. No.: #131	HRP-Protector™ Cat. No.: #220
Blocking of non-specific binding sites on <i>surfaces</i> (e.g. plates)	+++	-	-	-	++	-	-
Blocking of HAMA and Heterophilic Antibody interferences, cross-reactivities <i>in solution/matrix effects</i>	-	+++	+++	+++	-	-	-
Stabilization of <i>immobilized</i> proteins / antibodies	+	-	-	-	+++	-	-
Stabilization of proteins / antibodies <i>in solution</i>	-	+	++	+++	-	+++	++
Preserving of antibody-HRP / fluorophore conjugates	-	-	+++	+++	-	++	+++